

MECTA ELECTRONIC MEDICAL RECORD

MECTA EMR^{©2007} FOR SPECTRUM ECT



Prepare For Paperless Medicine From MECTA

The MECTA EMR[©] with its elegant simplicity of use is an extraordinarily powerful software program allowing extensive organization and analysis of your ECT treatment data.

Government and hospital organizational mandates require digitization of patient medical records. Improvement of patient safety and services, patient privacy, and cost and time efficiencies are ensured with MECTA's newest software program.

Automatically records up to 52 data points, six channels of physiological monitoring, and allows manual inputting of 700 data points.

All physiological monitoring and data can be printed and exported to text or Excel files for further analysis.

Network with hospital Electronic Health Record database.

Ensure time and cost efficiencies with easy data entry that includes new drop down menus, pop up menus, and autofill features.

Merge, Import, Export, Query, Backup and Restore databases.

AUTOMATICALLY GENERATED REPORTS

ECT Session Note Report

ECT Session Report

ECT Stimulus and Monitoring Report

History of ECT Treatments

MANUAL REPORTS


ECT Evaluation

Pre-Course Anesthesia Evaluation

Pre-ECT Nursing Report

ECT Anesthesia Report

Post-ECT Nursing Report



MECTA's Newest Software Product is NOW AVAILABLE!

■ This Sequel Server[™]-based, networkable MECTA EMR[©] will revolutionize the data acquisition of your ECT treatments. The SPECTRUM alphanumeric data, physiological monitoring and up to 700 data points automatically and manually entered by psychiatrists, anesthesiologists and nurses complete ONE ECT Session Note Report and/or FIVE Manual Reports. FOUR automatically generated Reports provide data for review. All can be saved as pdfs on your pc, laptop, hospital server or hospital Electronic Health Record.

Current requirements:

Recommended Computer System: Pentium processor @ 600 MHz, 2 GB RAM, and 100 GB hard drive space. DVD RW drive. External backup to a separate hard drive or compact disk. System software: Windows XP Pro SP3 or Windows 7 PRO.



Visit our website at www.mectacorp.com

Please contact MECTA representative in your area for information and upgrades.